


## TECHNICAL SHEET

<b>Reference</b>	<b>IN01/02</b>
<b>Trade Name:</b>	<b>pH Indicator strips</b>
<b>Product Description</b>	pH Indicator strips for human gastric aspirate testing
<b>Intended Use</b>	The pH Indicator strips are used to determine the correct placement/misplacement of nasogastric feeding tubes by testing the pH value of human aspirate.
<b>Raw Materials:</b>	1. Pure paper ( Cellulose) 2. Indicator dye



<b>Packaging</b>	Multiple package (see table below). QML (minimum delivery quantity): 1 Box  <b>Reference Features: Unit</b> <table border="1" style="width: 100%;"> <thead> <tr> <th>code</th> <th>Box/Tube</th> </tr> </thead> <tbody> <tr> <td>IN01</td> <td>200 per box</td> </tr> <tr> <td>IN02</td> <td>100 per tube</td> </tr> </tbody> </table>	code	Box/Tube	IN01	200 per box	IN02	100 per tube
code	Box/Tube						
IN01	200 per box						
IN02	100 per tube						
<b>Shelf Life:</b>	Non sterile						
<b>Usage</b>	Single use device						
<b>Storage</b>	Normal storage procedures and conditions						
<b>Disposal</b>	The user must follow the legal regulations regarding disposal of hospital waste.						
<b>Manufacturer</b>	Medicina ltd Units 1-4 Rivington View Business park Station Road Blackrod						

	Bolton BL6 5BN
<b>Classification</b>	Class I
<b>Conformity Assessment Route:</b>	MDD 93/42/EEC, Annex VII